

**Neuroblastoma dissociation technique**

Date: \_\_\_\_\_ Tumor Name: \_\_\_\_\_

Diagnosis:

Clinical History:

Procedure:

Surgeon:

Pathologist:

Neuroblastoma Tumor Processing – Cell Dissociation with RBC lysis:

1. Tumor rinsed with PBS<sup>-/-</sup>; use sterile scalpels to chop/mince tumor into small pieces/slurry
2. Transfer to 50ml conical tube (s) → use multiple tubes as needed for adequate dissociation if tumor > 5 grams or dissociate in 200 ml flask.
3. Fill with PBS<sup>-/-</sup> to 50ml mark
4. Dissociation
  - a. Remove Trypsin from -20 freezer and thaw in 37 degree water bath
  - b. Add trypsin to each 50 ml conical tube or flask with PBS and tumor
    - i. Tumor <5g = 600uL; Tumor 5-10g = 900uL; Tumor >10g = 1200uL
  - c. Place in warm 37 degree water bath x 10 minutes
  - d. Check to see if dissociated properly (goal is “slimy” when pipetted near tumor specimen). If too thick add more PBS<sup>-/-</sup> and incubate again.
    - i. Total time in water bath \_\_\_\_\_
  - e. Once complete stop dissociation by adding to the tube with the tumor i. STI (equal amount to trypsin added above)
    1. Tumor <5g = 600uL; Tumor 5-10g = 900uL; Tumor >10g = 1200uL
    - ii. DNase – 1/10<sup>th</sup> volume of STI/Trypsin
      1. Tumor < 5g = 60uL; tumor 5-10g = 90uL; tumor > 10g = 120uL
    - iii. 1M Magnesium Chloride – 1/10<sup>th</sup> volume of STI/Trypsin
      1. Tumor < 5g = 60uL; tumor 5-10g = 90uL; tumor > 10g = 120uL
    - iv. Add DNase and MagCl in 30uL increments (each) until “slime” (DNA) resolves; mix well and allow to sit at room temp 2-3 min between each addition.
5. Filter with a 40uM cell strainer (with lip) over the erylenmeyer flask = \_\_\_\_\_ mL cell suspension
6. View under microscope to determine viability
  - a. If many clumps, use 20ml syringe and 22 gauge needle to draw up suspension and transfer to another flask
7. Spin cells at 450g (G=RCF) x 5 minutes, aspirate and discard supernatant
8. Add 10ml of RBC lysis solution to each tube with cells.
9. Incubate at room temp x 10 minutes.
10. Mix solution of PBS/10% FBS: add to each tube to a total of 50ml.
11. Spin at 450g x 5 minutes, aspirate and discard supernatant.
12. Resuspend cell pellet in PBS/10%FBS.
13. Count cells – if too many large clumps, use 20guage syringe and 18 guage needle

